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009575EX-20191018

bre.flex 2nd Edition für thermopress 400

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(EN) Instructions for use

bre.flex 2nd Edition for thermopress 400

Before using this product, please read the instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

bre.flex 2nd Edition is a thermoplastic resin based on polyamide (PA) which can be processed with the „thermopress 400“ injection moulding system of bredent Co. The processing parameters to ensure reproducible quality are preset in the thermopress 400 injection moulding system.

bre.flex 2nd Edition is a biocompatible resin that is used to fabricated metal-free, removable restorations and splints.

For the approved dental indications, please refer to the table listed under item 2. bre.flex 2nd Edition is available in 3 different colors (clear, PC20, pink veined) and is prefilled in aluminium cartridges as a granular material available in different filling quantities or in a bulk quantity of 500 g.

2. Indication

Product	Removable/ partially removable	fixed	Additional information
bre.flex 2nd Edition (polyamide) semi-flexible	Clasp, attachment and telescopic injection-type bases, full and partial dentures, splints, sports mouthguards	not recommended	Colors: clear, PC20, pink veined

2.1 Target patient group

The patient group consists of adults with residual teeth in the upper and lower jaw.

2.2 Intended users

The product may only be used by dentists, dental technicians and adequately trained dental staff.

3. Contraindication

Crowns and bridges

The product must not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

Material samples for allergy sufferers:

Material sample	REF	Color	Qty	Type
bre.flex 2nd Edition	99253EF8	clear	1 piece	Polyamide
bre.flex 2nd Edition	99253EF7	PC20	1 piece	Polyamide
bre.flex 2nd Edition	99253EF6	pink veined	1 piece	Polyamide

Thermopress materials – allergenic substances

The following information is not binding and based only on professional assessment. It is not based on results of analyses or information provided by our raw material manufacturers or suppliers. The respective material should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

bre.flex 2nd Edition (polyamide)		
can be contained in traces:	yes	no
Methyl methacrylate		x
Ethylene glycol dimethacrylate (EDGMA)		x
Triethylene glycol dimethacrylate (TEGDMA)		x
Urethane dimethacrylate (UDMA)		x
Bis-GMA (Bisphenol A-dimethacrylate)		x
2,2-Bis-(2-hydroxy-3-methacryloxy-propoxy)-phenyl]propane		x
Butane diol dimethacrylate (BDDMA)		x
Camphor quinone		x
Benzoyl peroxide (BPO)		x
N,N-dimethyl-p-toluidine		x
N,N-Bis-(2-hydroxyethyl)-p-toluidine		x
Hydroxy methoxy benzophenone		x
Hydroquinone		x
Hydroquinone monomethyl ether (HQME)		x
Butyl-hydroxy-toluene (BHT, DBPC)		x
Formaldehyde		x
Colophony		x
Cobalt chloride		x
Phthalates		x

4. Hazard and safety information

Do not inhale vapors during the melting phase. Use caution in case of spilled granular material! Risk of accident by slipping! For further safety information about the product, please refer to the safety datasheet. The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

5. Storage and durability information

The granular material needs to be stored under clean and dry conditions. If stored correctly, the material has an unlimited shelf life. Due to the moisture sensitivity of the granules, each cartridge is individually sealed into an aluminium packaging. Only sealed cartridges may be used! Cartridges without packaging cannot be stored and must be disposed of immediately.

Information - 500 g bulk quantity:

Please ensure dry storage to prevent the material from absorbing any moisture. Water vapor resulting from residual moisture during the heating phase of the cartridge in the unit may result in bursting of the cartridge! Prior to the use, please predry the non-mixed granular material at least 2 hours at a temperature of 80°C. To maintain consistent quality of the granular material, the integration of foreign substances or contaminated particles must be avoided when filling the cartridges. Different thermoplastic granular materials must not be mixed with one another!

6. Processing/handling

Model fabrication:

When processing bre.flex 2nd Edition, a formaldehyde-free class IV plaster is used for the fabrication of plaster models. Exakto-Rock 5 (REF 5700SE50 ivory or 5700S850 white) (see table) is particularly suitable.

Mixing ratio:

Use of formaldehyde-free class IV plaster

Material	Model	Exakto-Rock 5	distilled water	Hardening time
bre.flex 2nd Edition	per	100 g	20 ml	1 h

Please invest the finished model into the flask once the total expansion time for the plaster has elapsed.

Wax model:


The dimensions used for wax models should be very generous. In other words, clasps and other delicate elements should be made thicker. A modelling wax with a low melting point, e.g. Biotec modelling wax (REF 51000610), is ideal for the boiling out phase with boiling water. As a result, any color residues on the injection-moulded thermoplastic part are avoided.

Mechanical retentions:

To ensure reliable bonding between denture tooth and denture base, all contact surfaces of the prefabricated teeth need to be sandblasted with 110 µm aluminium oxide.

To achieve reliable bonding of denture tooth and denture base, the contact surfaces of all prefabricated teeth need to be sandblasted with 110 µm aluminium oxide. A circumferential groove (1) needs to be milled in the cervical area and mechanical retention (undercut) (2) needs to be included in the basal area (see drawing).


The sandblasted teeth with mechanical retentions are additionally conditined with visio.link primer (REF VLPMMMA10) (observe instructions for use of visio.link).

Rebasing in the injection technique (retentive) is used for high-quality extensions or repairs.

Film sprue:

A film sprue must be attached to the wax model in order to benefit from the outstanding mechanical properties of thermoplastics. Apply this film sprue using a 1.5 mm thick wax plate. Using a 10 mm thick injection moulding channel (REF 43001720), attach the molten thermoplastic directly to the model via the 1.5 mm wax plate.

Injection parameters:								
Invest-ment technique	Material	Target temp. in °C	Heat-ing time in min.	Pressure time in sec.	Speed	Force	Flask temp. in °C	Program number in thermopress 400 Version 6.62
Flask	bre.flex 2nd Edition	280 °C	15 min.	90 sec.	6	165	60 °C	4

The following settings are required for compressed air-operated injection moulding systems:

Material	Parameter - pressure (force) in bar
bre.flex 2nd Edition	9,5

To facilitate opening the halves of the flask and to avoid damaging the moulded object, the press-out device and punch (REF 14000904) should be used. Any larger pieces of plaster must be carefully removed from the moulded object using a pneumatic vesting chisel. Fine cross-cut silicone burs (REF 5263QG60) have proven to be perfectly suitable for finishing and reducing the moulded object.

The special bur set (REF 33000830) is recommended for easier finishing of the thermoplastics.

Polishing

First of all, smooth the structure down using sandpaper with a 280 grain, before evenly smoothing it down with the handpiece using a 400 grain. Then the structure is prepolished with pumice powder and a goat-hair brush (REF 35000610) using the polishing motor without applying any pressure. Complete high gloss polishing with the handpiece using a goat-hair brush (REF 35000550) and the high-gloss polishing paste Abraso-Star K50 (REF 52000161) or Abraso-Starglanz (REF 52000163).

7. Technical data	
Composition	
Polyamide	> 99,9 %
Pigments	< 0,1 %
Properties according to DIN EN ISO 20795-1 (typical values)	
Flexural strength	65 MPa
Modulus of elasticity	1400 MPa
Water absorption	30 µg/mm³
Solubility	0 µg/mm³

8. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff. These instructions for use need to be handed over to the user as an accompanying document for insertion or further processing of the product. Only original tools and components may be used for processing.

CAUTION

Proper fit of the prosthetic restoration or the retention element must be checked when inserting the restoration and during the patient's annual recall appointment.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

The user and/or patient is obliged to report any serious incidents that have occurred in relation to the device to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is registered.

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