ZL Microdent Attachments

DE

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
ZL Microdent Attachments

EN

Instructions for the use of attachments, latch attachments, anchors, and retention elements. Prior to use, carefully read the instructions for use and instruction manual for individual products. For use only by a dental professional within the recommended indications.

Product Description

The ZL attachments are a precision mechanical construction element for incorporation into prostheses and serve to fix and stabilize removable dentures. They contribute to the functional restoration of the craniomandibular chewing system.

Attachment Variants	Product name	Application area
Intracoral Attachment	ANTERIOLOCK DUOLOCK ^{30°} DUOLOCK ^{90°} DUOLOCK ^{Prospektiv} DUOLOCK ^{75mm} INTRA LOCK LOGA MEGALOCK VENTRALOCK	Rigid support elements used in free-end and switching prostheses as well as removable bridges in patients with no or little atrophied alveolar ridge. The pull-off forces can be adjusted.
Extracoral Attachment	ACRYLOCK CADLOCK CENTRALOCK II	
Latch Attachment	CENTRALOCK MI ROBOLOCK	Rigid, extracoral latch attachment used in uni- and bilateral free-end and switching prostheses, removable bridges and hybrid prostheses.
Anchor	ANCHOR ATTACHMENT COMBIOLOCK ZL ANCHOR SYSTEM	Enables constructions for pro-thetic restoration with bilateral free-end and switching prostheses and removable bridges. Due to the individual bar length design, the combi bar is ideal for supplying variable switching gaps.
	WUKA	For the retention and stabilization of mucosal-supported overdenture prostheses with a small residual pillar stock (minimum number 2).
Ball Attachment	SPHÄROLOCK	
Retention element	SECURALOCK MULTIOCK	Retention element for double crowns with friction loss

Intended use

ZL attachments are encapsulated in dental restorations to anchor removable prosthesis in oral cavity.

Special Indications

- DUOLOCK^{75mm} is indicated only for use in highly sensitized patients.
- DUOLOCK^{Prospektiv} - The special construction of this attachment makes it possible that in the event of a possible foreseeable loss of pillars, the intracoral primary parts can be used as a base and anchoring element for fixed combined and removable dentures.
- VENTRALOCK - An extremely small precision attachment which may only be used for the indication in the anterior tooth area.
- CADLOCK - It is specifically used for CAD/CAM applications and for frameworks made of cobalt-chrome, titanium and zirconium.

Target patient population

The intended patient group are adults with remaining teeth in the upper/lower jaw and/or upper/lower jaws treated with implants.

Contraindications

- Patients with an existing allergy to one or more elements of the attachment materials.
- Inability or unwillingness of the patient to follow the aftercare instructions.
- Patients with extreme bruxism or further uncontrolled para-functional habits.

Cautions

- ZL Attachments may only be used by dentists and dental technicians who have appropriate experience, knowledge and training.
- Refer to the product-specific instruction manual for detailed information. The instructions in the corresponding manuals must be followed.
- The components of the ZL Attachments cannot be combined in conjunction with retaining elements from other manufacturers. Only original ZL instruments, tools and auxiliary parts must be used.
- The devices are delivered non-sterile.
- The products labelled with “Do not re-use” are intended for single use only. If the product is used more than once, the following risks may arise such as cross-contamination, malfunction, poor fit etc.
- At Occurrence of sensitivity reactions should discontinue use and consult a doctor.
- After incorporating the prosthetic restoration, regular control examinations are important as recommended by the physician.
- The device has not been evaluated for safety in the MR environment. Caution for MRI environment – to remove the removable part of the prosthesis before MRI scanning. When scheduling an MRI or other radiology services while wearing ZL attachments, the patient must always inform MRI or other radiology staff prior to the procedure so that proper coordination of care can be arranged.
- The components are delivered in a non-sterile state. The dentures must be thoroughly cleaned before they are placed in the patient’s mouth. Please follow the instructions for use of the manufacturer of the disinfectants.
- To ensure a secure fit of the prosthesis and to avoid overloading the attachment parts, the denture bearing must be checked at the six-monthly recall and the prosthesis must be relined if necessary.
- A prosthesis with loose, damaged or broken attachments should not be worn. Otherwise this may lead to choking, swallowing, aspiration or other injuries.
- The patient must be instructed on the correct insertion and removal of the prosthesis as well as daily prosthesis care.

Storage

The device must be stored in a dry and dust-free environment in its original packaging.

Disposal Instructions

As per applicable regional/national and international guidelines on dental waste management.

Reporting of undesired effects

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilanceZL@gc.dental. In this way you will contribute to improve the safety of this product.

Explanation of the symbols

	Manufacturer		Batch code
	Date of manufacture		Catalogue number
	Medical device		Keep dry
	Do not re-use		Information on packaging size
	Caution, consult accompanying documents		Consult the instructions for use

		
ZL MICRODENT Attachment GmbH & Co. KG - a GC Europe AG company - Harkortstraße 2 58339 Breckerfeld/Germany		

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